



FUNDAMENTALS OF LAPAROSCOPIC SURGERY™ TRAINER SYSTEMS



LEARN MORE AT

fls-products.com

PACKAGE SUPPLIED	3
TEST CERTIFICATION	4
SAFETY & CLEANING INSTRUCTIONS	6
SET UP	7
ATTACHING THE TROCARS	7
ATTACHING THE FRAME	8
ASSEMBLING THE TRAINER BOX	8
CONNECTING A MONITOR	9
TIPS	10
• ADJUSTING THE CAMERA TO CENTER THE IMAGE	10
REMOVE PLASTIC CAP	11

THREE SYSTEMS AVAILABLE





	CE	randi fazzle dal filmate
EC DECL	ARATION OF CONFORMITY	DECLA
	accordance with ISO IEC 17050-1 : 2010	In ac
	We at	
	Limbs & Things Ltd	
	Linus & Things Ltu	
	of	
Sus	ssex Street, St. Philips, Bristol, BS2 0RA	Su
Declare under	our sole responsibility that the product described below:	Declare unde
Brand and Model	FLS (Fundamentals of Laparoscopic Surgery) Trainer, Part	Brand and Model Nat
Name	numbers: 50301,50302, 50306	
	tested and is in conformity with the following standards:	Has been t
SUBJECT	STANDARD LVD Directive – 2014/35/EU	SUBJECT
Electrical Safety	EN 61010: 2010 + A1:2019 Clause	Electrical Safety
	5,6,7,8,9,10,11,12,13,14,15,16 & 17	
Electromagnetic	EMC Directive – 2014/30/EU	Electromagnetic
Compatibility	Conducted emissions:	Compatibility
	EN 55011:2009 + A1:2010	· · · · · · · · · · · · · · · · · · ·
	Radiated emissions:	
	EN 55011:2009 + A1:2010	
	Voltage Fluctuations (Composite Video and USB Mode):	
	EN 61000-3-3:2013	
	Electrostatic Discharge (Composite Video and USB Mode)	
	EN 61000-4-2:2009	
	EFT/Bursts (Composite Video and USB Mode)	
	EN 61000-4-4:2004 + A1:2010	
	Surge Immunity (Composite Video and USB Mode)	
	EN 61000-4-5:2006	
	Conducted Immunity (Composite Video and USB Mode)	
	EN 61000-4-6:2009	
	Dips & Interruptions (Composite Video and USB Mode)	
Podio porformanco	Dips & Interruptions (Composite Video and USB Mode) EN 61000-4-11:2004	Dadia auformana.
Radio performance RoHS	Dips & Interruptions (Composite Video and USB Mode)	Radio performance RoHS

UK CA RATION OF CONFORMITY ccordance with UK Government Guidance We at Limbs & Things Ltd of ussex Street, St. Philips, Bristol, BS2 0RA ler our sole responsibility that the product described below: ame FLS (Fundamentals of Laparoscopic Surgery) Trainer, Part numbers: 50301,50302, 50306 tested and is in conformity with the following standards: STANDARD Electrical Equipment (Safety) Regulations 2016 BS EN 61010: 2010 + A1:2019 Clause 5.6.7.8,9,10,11,12,13,14,15,16 & 17 Electromagnetic Compatibility Regulations 2016 Conducted emissions: BS EN 55011:2009 + A1:2010 Radiated emissions: BS EN 55011:2009 + A1:2010 Voltage Fluctuations (Composite Video and USB Mode): BS EN 61000-3-3:2013 Electrostatic Discharge (Composite Video and USB Mode) BS EN 61000-4-2:2009 EFT/Bursts (Composite Video and USB Mode) BS EN 61000-4-4:2004 + A1:2010 Surge Immunity (Composite Video and USB Mode) BS EN 61000-4-5:2006 Conducted Immunity (Composite Video and USB Mode) BS EN 61000-4-6:2009 Dips & Interruptions (Composite Video and USB Mode) BS EN 61000-4-11:2004 N/A The Restriction of the use of certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012 e equipment named above has been designed to comply with the relevant Signature eferenced Standards. lull ging Director & Things, Bristol, England 05/11/21 Date

FC <u>DECLARATION OF CONFORMITY</u>			
Sussex	Street, St. Philips, Bristol, BS2 0RA		
Declare under our s	ole responsibility that the product described below:		
Brand and Model Name	FLS (Fundamentals of Laparoscopic Surgery) Trainer, Part numbers: 50301,50302, 50306		
Has been tester	d and is in conformity with the following standards:		
Technical regulations	FCC CFR47 Part 15:107 & 15: 109 Class B limits		
Test laboratory	Kiwa Electrical Compliance Unit 8 Woodfieldside Business Park Pontllanfraith Blackwood United Kingdom NP12 2DG		
Test report number	R20-5155		
Test report date	07th April 2020		
Signed Nick Hull Position Managing Dir Location Limbs & Thin	Signature ector gs, Bristol, England Date 13 November 2004		

SAFETY & CLEANING INSTRUCTIONS

READ THE SAFETY INSTRUCTIONS CAREFULLY BEFORE USING YOUR FLS™ TRAINER SYSTEMS





Follow all instructions. Keep these instructions. Heed all warnings.



imbs&

things

Do not install near any heat sources, do not expose the product to flames or use near naked flame sources. Do not use flammable solvents near or on this product.

Do not use any attachments

other than those supplied by

_ _ _ _ _ _ _ _ _ _ _ _ _ _ _ _

the manufacturer.

DON'Ts

The product must be used by appropriately trained adults only.



Do not use the product near water, this product contains electronics.





Only use this product indoors.

Clean the product with a damp soft cloth or sponge, using only warm water with mild detergent.



Do not disassemble the product or attempt to repair a faulty unit, please return the product to Limbs & Things for assessment and repair.

_ _ _ _ _ _ _ _ _ _ _ _ _ _ _ _ _





If a product is damaged or worn, it should be repaired by Limbs & Things. To arrange this, please contact your point of purchase who will be able to provide a repair quote.



The product is not a toy and must not be used by children or untrained staff.

STORAGE & TRANSPORTATION



Always ensure that the products are properly packed and secured during transportation and storage in order to prevent personal injury or damage to the product.



ATTACHING THE FRAME





ASSEMBLING THE TRAINER BOX





CONNECTING THE MONITOR



Turn TV on ·····

If monitor is not displaying the correct channel

2 Press 'Source' and select AV





After a few seconds the AV channel will display the camera view TIPS

ADJUSTING THE CAMERA TO CENTER THE IMAGE





REMOVE PLASTIC CAP

Ensure the cap covering the camera lens is removed before use

















Part No: 065-701

Issue 13 - November 2021 © 2021 Limbs & Things