

Product Recycling Information: CMPT Product Range

Below is a pictorial guide to components within the Limbs & Things CMPT product range, their material/s and recyclable status.

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Part Description	Model	Material	Recyclable	Disposal Information
Carry Cases 	All CMPT Product No. 40209/40210	PET, Foam		Widely recycled
CMPT Modules 	All CMPT Modules	Mixed (Silicone, Latex)		Specialist recycling or landfill
Lubricant Packets 	All CMPT	Mixed		Specialist recycling or landfill
		Lubricant (water soluble)		Lubricant - safe to dispose
CMPT Foreskins 	All CMPT Product No. 60961/60974	Silicone		Specialist recycling or landfill
CMPT Pubic Bone 	All CMPT Product No. 60960	ABS		Recyclable, non-standard recycling

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

CMPT Lower Torso 	All CMPT	Silicone		Specialist recycling or landfill
CMPT Abdominal Insert 	All CMPT Product No. 60959/60973	Silicone, Foam		Specialist recycling or landfill
CMPT Base 	All CMPT	Silicone		Specialist recycling or landfill
Packaging 	All CMPT	Cardboard		Widely recycled
LED Pen Torch 	Adv CMPT 60951/60965	Mixed		Specialist recycling - e-waste & battery recycling facilities
CMPT Augmented Reality Mat 	Adv CMPT Product No. 60135	Digiflor Tactile fire-rated slip-rated cushioned vinyl flooring 1000gsm		Recyclable, non-standard recycling

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.