

Surgical Female Pelvic Trainer (**SFPT**) Mk 2


60251

The Surgical Female Pelvic Trainer (**SFPT**) is a simulator for training in a wide range of laparoscopic gynecological surgical techniques and procedures. The uterus component can be easily removed and replaced after use.



Compatible for use
with harmonic scalpels

Realistic bleeding
when arteries are cut




Learning has never been more lifelike

limbsandthings.com



Abdominal Skin
with realistic tissue
properties for
trochar insertion



Working with gynecology
instruments

Product Qualities

Anatomy

- Includes:
 - Uterus
 - Cervix
 - Bladder
 - Ovaries
 - Ureters
 - Ovarian artery
 - Dunn's pouch
 - Ovarian cyst
 - Fibroid
 - Cardinal ligament
 - Vagina
 - External introitus
 - Fallopian tube
 - Broad ligament
 - Uterine artery
 - Pouch of Douglas
 - Bowel and sigmoid
 - Ectopic pregnancy
 - Uterosacral ligament
- Abdominal skin with realistic tissue properties for trochar insertion
 - Pelvic floor with sympathetic features

Realism

- Insufflation of abdomen is represented

Key Features

- Provides practice and role playing for gynecological perioperative teams
- Includes a realistic vaginal canal and cervix for inspection, in which a uterine manipulator can be used
- Realistic bleeding when arteries are cut

Versatility

- Compatible for use with harmonic scalpels

Safety

- Capsule in uterus contains latex

Skills

- All skills relevant to performing the following procedures:
 - Salpingectomy
 - Myomectomy
 - Hysterectomy
- Salpingostomy
- Cystectomy
- Oophorectomy
- Dissection down to and location of ureters
- Insertion and use of uterine manipulator
- Working with various gynecology instruments
- Identification of anatomical landmarks

- SFPT Mk 2 Abdomen
- 60252 SFPT Mk 2 Surgical Uterus
- 60253 SFPT Mk 2 Abdominal Skin Pad
- 60254 SFPT Mk 2 Perineum
- 60255 SFPT Mk 2 Bowel
- 00020 Venous Blood Starter Kit
- SFPT Mock Blood Kit (includes reservoir and pump)
- Blunt dispensing needle and syringe
- Bed securing strap

Consumables

- 60256 SFPT Uterus (x5)
- 60257 SFPT Uterus (x10)

Curriculum Reference

ACGME OB/GYN Milestones (p. 12):

Gynecology Technical Skills: Laparotomy (e.g., Hysterectomy, Myomectomy, Adnexectomy): Applies innovative and complex approaches to laparotomy and implements treatment plans based on emerging evidence.



For demonstration or to order

Limbs & Things Inc. P.O. Box 15669, Savannah, GA 31416, USA

E: customerservice@limbsandthings.com T: (912) 629-0357

limbsandthings.com