



MALE RECTAL EXAMINATION TRAINER

STANDARD & ADVANCED

This brand new trainer offers realistic, repeatable training in the core skill of digital rectal or PR examination.

Both models offer a contractible anal sphincter, allowing for anal tone assessment and a faecal impaction clip in module.

60170	60186	STANDARD
60171	61024	ADVANCED



Learning has never been more lifelike

limbsandthings.com


The **Standard** model focuses on learning core rectal nursing procedures and excludes the pathologies.



The **Advanced** model comes with a range of prostates and 2 perineums, allowing trainees to become familiar with normal and abnormal findings and acquire key diagnostic skills.

External ano-rectal illustrations are ideal for incorporation into scenarios and OSCE assessments.

PATHOLOGICAL PERINEUM




NORMAL

UNILATERAL ENLARGED


BILATERAL ENLARGED


UNILATERAL CARCINOMA

BILATERAL CARCINOMA



EXTERNAL CONDITIONS ILLUSTRATION PACK





ANAL TONE ASSESSMENT

QUALITIES

* Advanced model

ANATOMY

- Buttocks, anus, rectum, prostate and perineum
- Pathological perineum – with polyp and rectal cancer*
- 4 abnormal prostates*:
 - Benign unilateral enlarged
 - Benign bilateral enlarged
 - Unilateral carcinoma
 - Bilateral carcinoma

REALISM

- Soft, partable buttocks
- Realistic anus with resting tone and ability to contract

KEY FEATURES

- Simulation of sphincter contraction allows trainees to assess anal tone
- Interchangeable prostates can be quickly and easily inserted, out of sight of the trainee*
- Left lateral positioning
- Addition of impacted faecal matter allows trainees to recognise and distinguish this common finding
- High quality illustration pack of 9 external anal conditions*

VERSATILITY

- Trainer can also be presented in 'semi-standing' position using the optional **Standing Position Stand (60172)**

SAFETY

- Latex free

OPTIONAL EXTRAS

- **60172** Standing Position Stand
- **40209** Medium Carry Case

SKILLS

- Digital examination of the anus, rectum and prostate
- Assessment of anal tone
- Identification of faecal matter in the rectum
- Identification of external ano-rectal conditions*
- Professional to patient communication

CONTAINS

	60170	60171
Base Unit	✓	✓
Left Lateral Stand	✓	✓
Normal Perineum with Anal Tone	✓	✓
Pathological Perineum Anal Tone	✗	✓
Normal Prostate	✓	✓
Benign Unilateral Enlarged Prostate	✗	✓
Benign Bilateral Enlarged Prostate	✗	✓
Unilateral Carcinoma Prostate	✗	✓
Bilateral Carcinoma Prostate	✗	✓
Impacted Faecal Matter	✓	✓
External Conditions Illustration Pack	✗	✓
Lubricant	✓	✓





CURRICULUM REFERENCES:

UNDERGRADUATE MEDICAL

Medical Deans' Clinical Assessment Blueprints for the medical graduate

4. Conducts ... organ(s)/system specific physical examination as appropriate. 5. Integrates the information obtained to arrive at an appropriate diagnosis and differential diagnosis.

POSTGRADUATE MEDICAL – JUNIOR DOCTOR

CPMEC Australian Curriculum Framework for Junior Doctors v 3.1, 2012, p.7

History and Examination: Performs a comprehensive examination of all systems; Elicits symptoms & signs relevant to the presenting problem or condition. Discriminates between the possible differential diagnoses relevant to a patient's presenting problems or conditions

GP TRAINING

RACGP Curriculum for Australian General Practice 2016, CS16 Core Skills unit, p.23

CS2.2.2.2a Demonstrate appropriate and respectful physical examination technique

CS2.2.2.5a Identify situations where physical examination findings are inconsistent with history

CS2.2.2.2b Effectively utilise appropriate clinical tools to optimise examination

CS2.2.2.3b Effectively summarise key examination findings for the patient



FOR DEMONSTRATION OR TO ORDER

Limbs & Things Australia Pty Ltd. 92 Wedgewood Road, Hallam, Victoria 3803, Australia

E: info@limbsandthings.com T: 00 61 3 9708 6511

limbsandthings.com