




# Pediatric FAST/Acute Abdomen Phantom

## Ultrasound training for pediatric FAST procedures

This world's first pediatric ultrasound torso phantom provides opportunities of hands-on training in ultrasound examination of 2-year-old patient

### Product Supervision

Shunsuke Nosaka, MD, PhD  
 Director of Radiology  
 National Medical Center for Children and Mothers




### FEATURES

- 1 | The phantom of a 2-year-old, thoracoabdominal organs, bone structure, free-fluid to learn FAST procedures and pathologies that are commonly seen in pediatrics.
- 2 | With this phantom, trainees can acquire skills in the basics of pediatric abdominal ultrasound as well as FAST/Acute sonography.

### SKILLS

- FAST procedures:  
Internal hemorrhage at perihepatic, perisplenic, pelvis, and pericardium area
- Sonography for acute pediatric patient:  
Bowel intussusception, appendicitis and biliary dilatation

### IMAGES



Pleural hemorrhage



Hydronephrosis



Appendicitis

### DESCRIPTIONS

#### SPECIFICATIONS

Size: W41 x D15 x H15 cm      Weight: 6.5kg/ 14.3lbs  
 W16 x D6 x H6 inch

#### SET INCLUDES

- 1 | ultrasound phantom
- 1 | storage case
- 1 | tutorial manual (DVD)

#### MATERIALS

Polyurethane elastomer  
 Latex free

#### REPLACEMENT PARTS

-

#### RECOMMENDED DEVICES

-

\* Ultrasound scanner is not included.

