

Product Recycling Information: Male Rectal Examination Trainer

Product Range

Below is a pictorial guide to components within the Limbs & Things Male Rectal Examination Trainer product range, their material/s and recyclable status.

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Part Description	Model	Material	Recyclable	Disposal Information
Comfigel 42g Tube Lubricant 	All Rectal	Tube (dependant on material)		Recyclable, check supplied material for details
		Lubricant (water soluble), cardboard		Lubricant - safe to dispose Cardboard - Widely recycled
Packaging 	All Rectal	Cardboard		Widely recycled
Plain Poly Bag 	All Rectal	Polyurethane		Widely recycled
Medium Carry Case 	All Rectal Product No: 40209	PET		Widely recycled

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Part Description	Model	Material	Recyclable	Disposal Information
MRE Trainer Base Unit 	All Rectal Product No. 60177/60191	Mixed		(See component breakdown below)
MRE Hips		Mixed		Special recycling or landfill
MRE Skin		Silicone		Special recycling or landfill
MRE Torso		Icorene 3560 Icbag 20 P1200		Recyclable, non-standard recycling
MRE Trainer Perineum 	Normal Product No. 60175/60189 Pathological Product No. 60176/60190	Mixed		(See component breakdown below)
MRE Perineum Assy		Mixed		Special recycling or landfill
In Line Hose Barb 3.2mm N230-1		Dupont Zytel 101F NC010		Special recycling or landfill
MRE Frame		Sabir C6200 PC/ABS		Recyclable, non-standard recycling

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Enema bulbs size 2, pastel blue	Soft non-phthalate, biobased Sarolit Plastisol 900CM		Special recycling or landfill
Silicone Tube	Silicone		Special recycling or landfill
MRE Perineum			
Barb Reducer 6.4 To 3.2, White N-Grs6x4	Nylon		Special recycling or landfill
Rectal perineum packing tray 0.5mm white PVC	PVC		Recyclable, non-standard recycling
Prostates & Faecal Matter Set 			
Male Rectal Product No. 60178	Mixed		(See component breakdown below)
Box & Insert	Cardboard		Widely recycled
Prostate - Normal	Silicone		Special recycling or landfill
Prostate – Benign Unilateral Enlarged			
Prostate – Benign Bilateral Enlarged			
Prostate – Unilateral Carcinoma			
Prostate – Bilateral Carcinoma			

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Impacted Faecal Matter		RTV Rubber A40		Special recycling or landfill	
					
Standing Position Stand		Product No. 60172	Mixed		(See component breakdown below)
					
3M Bumpon SJ5312		Rubber		Special recycling or landfill	
Stand		White Satin Acrylic Capped ABS		Recyclable, non-standard recycling	

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.