

Product Recycling Information: PROMPT Flex – Birthing Simulator

Product Range

Below is a pictorial guide to components within the Limbs & Things PROMPT Flex – Birthing Simulator product range, their material/s and recyclable status.

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Single Components:

Part Description	Model	Material	Recyclable	Disposal Information
Birthing Lubricant Gel Bottle 	All PROMPT Flex Product No. 10199	PP		Widely recycled
Bed Straps 	All PROMPT Flex Product No. 80130	Polyester		Widely recycled
Mini Carry Case 	All PROMPT Flex Product No. 40211	PET		Widely recycled
Large Carry Case 	All PROMPT Flex Product No. 01000			

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Packing Box & Cardboard Tube 	All PROMPT Flex	Cardboard		Widely recycled
Placenta - PROMPT Flex 	All PROMPT Flex Product No. 80123	Silicone		Specialist recycling or landfill
Abdomen 	All PROMPT Flex Product No. 80124/80224	Silicone		Specialist recycling or landfill
Perineum & Birth Canal 	All PROMPT Flex Product No. 80125/80225	Silicone		Specialist recycling or landfill
Cervix 	All PROMPT Flex Product No. 80126	Silicone		Specialist recycling or landfill
Pelvic Ring Clamp 	All PROMPT Flex Product No. 80107	Bayblend T65 XF PC/ABS		Recyclable, non-standard recycling
PROMPT Flex Augmented Reality Mat 	All PROMPT Flex Product No. 60133	Digiflor Tactile fire-rated slip-rated cushioned vinyl flooring 1000gsm		Recyclable, non-standard recycling

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

PROMPT Flex Baby - Standard 	PROMPT Flex Standard Product No. 80121/80221	Silicone		Specialist recycling or landfill
---	---	----------	---	----------------------------------

Composite Components:

Part Description	Model	Material	Recyclable	Disposal Information
Birthing Mother 	All PROMPT Flex Product No. 80120/80220	Mixed		(See component breakdown below)
PROMPT Flex Mother Upper & Lower		Mixed		Specialist recycling or landfill
PROMPT Flex APC & Handle		T65 Bayblend PC/ABS		Recyclable, non-standard recycling
PROMPT Flex Pelvic Ring - Toggle				
PROMPT Flex Hip Ball Socket		Acetal		Recyclable, non-standard recycling
25mm Webbing with D Ring		Polyester		Widely recycled
PROMPT Flex Pubic Bone		Silicone, Nylon, Zinc, and Stainless Steel		Specialist recycling or landfill

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

PROMPT Flex Pelvic Floor Plate	Mixed		Specialist recycling or landfill
M8 Star Knob	Zinc Plated Steel		Widely recycled
Washer-M4 Shakeproof	Stainless Steel		Widely recycled
Polyplast 30 Torx Pan 3x10mm			
Spring Plunger Delrin			
M4x20mm Penny Washers - A2			
Polyplast 30 Torx Pan 4x14mm			
Polyplast 30 Torx Pan 4x10mm			
Polyplast 30 Torx Pan 3x12mm			
M8x30mm Socket Button Screw			
Durable Dot Stud 550SST			
Durable Dot 7.5mm Eyelet Post (NPB) 550BCC0	Nickle Plated Brass		Widely recycled
A2 M6X12mm Button Head Screw (Bolt)	Steel		Widely recycled
M8x40mm Repair Washer	Aluminium		Widely recycled

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

PROMPT Flex IM Pad	PU Foam 73-78 SHORE 00		Recyclable, non- standard recycling
M8 x 18mm Nylon Washer	Nylon		Specialist recycling or landfill
Washer 4.3 X 9 X 1.4 Nylon			
PROMPT Flex APC Skin Peg	Silicone/Silicone Rubber		Specialist recycling or landfill
3M Bumpon SJ5312			
PROMPT Flex Pelvic Floor			
PROMPT Flex Upper Leg			
PROMPT Flex Leg End Cap			
Baby Outer			

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

e-waste Components

Please note: Due to the complexity of the PROMPT Flex Baby - Advanced, the unit is classified as e-waste. At the end of the units’ life, it should be disposed of at an appropriate e-waste facility, if one is available to you.

Below is a breakdown of the constituent parts for reference only. Do not dismantle this item, this action will void your product warranty and could cause injury to your person or property.

Part Description	Model	Material	Recyclable	Disposal Information
PROMPT Flex Baby – Advanced (Force Monitoring) 	PROMPT Flex Advanced Product No. 80122FM/80222FM	Mixed	 (See component breakdown below)	Specialist recycling - e-waste & battery recycling facilities)
 OTHER Sabic Lexan 940 PC, Polycarbonate	 Polyplast 30 Torx Pan 3x10mm, Stainless Steel	 Electronics, DC power cables	 Silicone/Silicone Rubber	

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.