

# Surgical Female Pelvic Trainer (SFPT)

Part No: 60251

**Learning Has Never  
Been More Lifelike**



Includes a realistic vaginal canal and cervix for inspection, in which a uterine manipulator can be used




Abdominal Skin with realistic pressure needed to insert Cannula

Insufflation of abdomen is represented



Compatible for use with harmonic scalpels



Realistic bleeding when arteries are cut



Provides practice and role playing for gynaecological perioperative teams

Anatomy includes:  
Uterus, Vagina, Cervix, External Introitus, Bladder, Fallopian Tube, Ovaries, Broad Ligament, Ureters, Uterine Artery, Ovarian Artery, Pouch of Douglas, Dunns Pouch, Bowel and Sigmoid, Ovarian Cyst Ectopic Pregnancy, Fibroid, Uterosacral Ligament and Cardinal Ligament

# Surgical Female Pelvic Trainer (SFPT)

Part No: 60251

The new Surgical Female Pelvic Trainer (SFPT) is a simulator for training in a wide range of laparoscopic gynaecological surgical techniques and procedures. The uterus component can be easily removed and replaced after use.

## Skills Gained:

- Salpingectomy
- Salpingostomy
- Myomectomy
- Cystectomy
- Hysterectomy
- Oophorectomy
- Dissection down to and location of ureters
- Insertion and use of uterine manipulator
- Working with various gynaecology instruments
- Identification of anatomical landmarks



- |   |                        |
|---|------------------------|
| - Uterus  | - Vagina               |
| - Cervix  | - External Introitus   |
| - Bladder   | - Fallopian Tube       |
| - Ovaries   | - Broad Ligament       |
| - Ureters   | - Uterine Artery       |
| - Ovarian Artery  | - Pouch of Douglas     |
| - Dunns Pouch   | - Bowel and Sigmoid    |
| - Ovarian Cyst  | - Ectopic Pregnancy    |
| - Fibroid   | - Uterosacral Ligament |
| - Cardinal Ligament   |                        |
| - Abdominal Skin with realistic pressure needed to insert Cannula |                        |
| - Pelvic Floor with sympathetic features                          |                        |

## Product Qualities



- Insufflation of abdomen is represented
- Provides practice and role playing for gynaecological perioperative teams
- Includes a realistic vaginal canal and cervix for inspection, in which a uterine manipulator can be used
- Compatible for use with harmonic scalpels
- Realistic bleeding when arteries are cut



## You'll Receive:

- 60252 1 SFPT Surgical Uterus
- 60253 1 SFPT Abdominal Skin Pad
- 60254 1 SFPT Perineum
- 1 Abdominal Support
- 1 Abdominal Shell

## + Consumables:

- 60256 SFPT Uterus (pack of 5)
- 60257 SFPT Uterus (pack of 10)



[limbsandthings.com](http://limbsandthings.com)

Limbs & Things Australia Pty Ltd, 92 Wedgewood Road, Hallam, Victoria 3803, Australia  
E: [info@limbsandthings.com.au](mailto:info@limbsandthings.com.au) T: +61 3 9708 6511 F: +61 3 9708 6566