





PARACENTESIS TRAINER

This trainer fulfills the fundamental component of core medical training in diagnostic and therapeutic techniques of Paracentesis, for either palpation or ultrasound guided techniques.

Trainees are able to identify the echogenic anatomy enabling them to learn how to safely insert a needle or catheter into the peritoneal cavity.

60100 LIGHT

60111 DARK



Learning has never been more lifelike

limbsandthings.com



QUALITIES

ANATOMY

- Torso featuring bony landmarks (pubis symphysis, iliac crest) and umbilicus
- Internal anatomy includes:
 - Liver
 - Bowel
 - Spleen
 - Floating Bowel

REALISM

- Realistic tissue and needle response

KEY FEATURES

- Palpation or ultrasound techniques can be practised (side by side)
- Internal echogenic anatomy to allow recognition of landmarks under ultrasound
- Two 3.5 litre chambers can be filled with water for practise drainage
- Cost effective consumable pads

VERSATILITY

- Self-sealing pads, withstanding up to 200 needle or up to 100 rocket catheter insertions
- Ability to insert and remove a drain
- Positions both supine and lying on side

SAFETY

- Latex free

CLEANING & CARE

- Skin surface is washable using soap and water

SKILLS

- Familiarity with the abdominal regions and underlying anatomy
- Palpation of anatomical landmarks
- Identification of excess fluid
- Through ultrasound, identification of safe needle/catheter insertion sites
- Insertion of needle into peritoneal cavity for therapeutic or diagnostic purposes
- Professional-to-patient communication

WHAT'S INCLUDED

60100

60111

Paracentesis Skin
 Paracentesis Pads (x2)
 Paracentesis Liver
 Paracentesis Spleen
 Paracentesis Bowel Insert
 Paracentesis Floating Bowel
 Paracentesis Clamping Frame
 Paracentesis Large Volume Refill System
 Base unit
 Carry case

CURRICULUM REFERENCES

EMERGENCY MEDICINE

Australasian College for Emergency Medicine Curriculum Framework, 2015, p. 74

JUNIOR DOCTORS

Confederation of Postgraduate Medical Education Councils, Australian Curriculum Framework for Junior Doctors, 2009

ONCOLOGY TRAINEES

The Royal Australasian College of Physicians: Medical Oncology Advanced Training Curriculum, 2013, p. 24

GENERAL SURGERY

General Surgery Curriculum, Royal Australasian College of Surgeons, General Surgeons Australia & New Zealand Association of general Surgeons, 2013

AUS V.3 Mar 2021



REQUEST A QUOTE



REQUEST A DEMO



BUY ONLINE



CUSTOMER SUPPORT

limbsandthings.com