

Product Recycling Information: Venipuncture Arms Product Range

Below is a pictorial guide to components within the Limbs & Things Venipuncture Arms product range, their material/s and recyclable status.

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Part Description	Model	Material	Recyclable	Disposal Information
Venipuncture Arm Skin 	All Arms Product No. 70311/70312	Silicone		Specialist recycling or landfill
ACF Vein Module 	Std & Adv Arms Product No. 70331/70332	Latex		Specialist recycling or landfill
Metacarpal Vein Module 	Std & Adv Arms Product No. 70333/70334	Latex		Specialist recycling or landfill
Basic Venipuncture Arm Vein 	Basic Arms Product No. 70335	Latex		Specialist recycling or landfill
Basic Venipuncture Arm Vein - Latex Free 	Basic Arms Product No. 70337	Silicone		Specialist recycling or landfill

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Lubricant Bottle 	All Arms	PP		Widely recycled
Carry Case 	All Arms Product No. 00301	PET		Widely recycled
Lube Sachet 	All Arms	Mixed		Specialist recycling or landfill
Concentrated Venous Blood Bottle 	All Arms Product No. 00020	HDPE		Widely recycled
Venipuncture Arm Infusion Tube 	All Arms Product No. 00335	Silicone/ Polyacetal resin		Specialist recycling or landfill
Venipuncture Arm Rest 	All Arms Product No. 00302	High Impact Polystyrene		Recyclable, non-standard recycling

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Part Description	Model	Material	Recyclable	Disposal Information
Pressurized Venipuncture Arm Mock Blood Supply 				
	Advanced Arms Product No. 70315	Mixed		(See component breakdown below)
EPDM Black Foam 1.5mm FT165E		Mixed		Widely recycled
Fluid Unit Case				
Bottle Plate - Lower		ABS		Recyclable, non-standard recycling
Fluid Unit Handle Tip Injection Moulded, Red				
Fluid Unit Valve Assy Complete		Nylon, PC Silicone PE		Specialist recycling or landfill
Bottle Assy		PE-HD, Silicone		Specialist recycling or landfill
Bottle Plate Assy		PP, Silicone		Specialist recycling or landfill
O Ring		Nitrile Butadiene Rubber (NBR)		Specialist recycling or landfill
1/8" Equal Female Connector		Aluminium		Widely recycled


Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Tee Tapper Connectors TES 4/R1/8	PBT, Brass		Specialist recycling or landfill
Hand Pump - Black V15600LBK	PVC		Recyclable, non- standard recycling
Panel Mount Coupler Nitto SPC-04PHB (NBR) Male Coupler	Polyacetal Resin		Specialist recycling or landfill
Panel Mount Coupler Nitto SPC-04PHB-BLU			
Silicone Tube	Silicone		Specialist recycling or landfill
3x10 Torx Pan SS 202101452K3010S	Stainless Steel		Widely recycled
Clear Feet 6mm Self Adhesive	Acrylic		Recyclable, non- standard recycling
Screw-3 X 16Mm-Pozi-Pan Head	Steel		Widely recycled
Nut-Full-Steel-M3-Z.P.			
0-10psi Gauge 0-10 PSI Panel Mount Gauge, 1/8 bsp, 40mm Diameter Black Steel Casing c/w Bracket			
Washer-M3-Int.Shake'Pf Zinc Plated	Zinc		Widely recycled
Black Foam Washer - 1.5mm	Foam		Widely recycled

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Refill Bottle				
	Advanced Arms Product No. 00303	Mixed		(See component breakdown below)
Bulkhead Nut		Stainless Steel		Widely recycled
Silicone Tube 4.8mmID x 8mmOD		Silicone		Specialist recycling or landfill
Bulkhead Connector BH4S4-4BN		Black Nylon		Specialist recycling or landfill
Quick Release Coupler Nitto SPC-06SH-IVR		Polyacetal resin		Specialist recycling or landfill
Quick Release Coupler Nitto (NBR) Female Coupler				
Venipuncture Arm				
	Standard & Advanced Arms Product No. 70316	Mixed		(See component breakdown below)
IV Arm Core		Expanded Polypropylene VB15 – 90G/L		Widely recycled
12mm Diameter Polypropylene Rod 100mm Length		PP		
IV Arm Skin Peg		Yungsox 5090T		
Silicone Tube		Silicone		Specialist recycling or landfill
IV Arm Sleeve				

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Medi Tube Clamp	PVC		Recyclable, non-standard recycling
IV Arm Basilic Vein	Latex		Specialist recycling or landfill
ACF Vein	Latex, Nylon		Specialist recycling or landfill
Gravity System Blood Supply 			
Standard Arms Product No. 70316	Mixed		(See component breakdown below)
Gravity Fed System Jug	Wonderlite PC-110		Widely recycled
Gravity Fed Fluid System Lid	TPE TC9GPN 89ShA		Recyclable, non-standard recycling
Gravity Fed System Stand	S235JR Mild Steel		Widely recycled

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.