



# MALE RECTAL EXAMINATION TRAINER

## STANDARD & ADVANCED

This brand new trainer offers realistic, repeatable training in the core skill of digital rectal or PR examination.

Both models offer a contractible anal sphincter, allowing for anal tone assessment and a faecal impaction clip in module.

60170	60186	STANDARD
60171	61024	ADVANCED



*Learning has never been more lifelike*

[limbsandthings.com](http://limbsandthings.com)



The **Standard** model focuses on learning core rectal nursing procedures and excludes the pathologies.




The **Advanced** model comes with a range of prostates and 2 perineums, allowing trainees to become familiar with normal and abnormal findings and acquire key diagnostic skills.


External ano-rectal illustrations are ideal for incorporation into scenarios and OSCE assessments.

PATHOLOGICAL PERINEUM



EXTERNAL CONDITIONS ILLUSTRATION PACK





ANAL TONE ASSESSMENT

QUALITIES

\* Advanced model

ANATOMY

- Buttocks, anus, rectum, prostate and perineum
- Pathological perineum – with polyp and rectal cancer\*
- 4 abnormal prostates\*:
  - Benign unilateral enlarged
  - Benign bilateral enlarged
  - Unilateral carcinoma
  - Bilateral carcinoma

REALISM

- Soft, partable buttocks
- Realistic anus with resting tone and ability to contract

KEY FEATURES

- Simulation of sphincter contraction allows trainees to assess anal tone
- Interchangeable prostates can be quickly and easily inserted, out of sight of the trainee\*
- Left lateral positioning
- Addition of impacted faecal matter allows trainees to recognise and distinguish this common finding
- High quality illustration pack of 9 external anal conditions\*

VERSATILITY

- Trainer can also be presented in 'semi-standing' position using the optional **Standing Position Stand (60172)**

SAFETY

- Latex free

OPTIONAL EXTRAS

- **60172** Standing Position Stand
- **40209** Medium Carry Case

SKILLS

- Digital examination of the anus, rectum and prostate
- Assessment of anal tone
- Identification of faecal matter in the rectum
- Identification of external ano-rectal conditions\*
- Professional to patient communication

CONTAINS

	60170	60171
Base Unit	✓	✓
Left Lateral Stand	✓	✓
Normal Perineum with Anal Tone	✓	✓
Pathological Perineum Anal Tone	✗	✓
Normal Prostate	✓	✓
Benign Unilateral Enlarged Prostate	✗	✓
Benign Bilateral Enlarged Prostate	✗	✓
Unilateral Carcinoma Prostate	✗	✓
Bilateral Carcinoma Prostate	✗	✓
Impacted Faecal Matter	✓	✓
External Conditions Illustration Pack	✗	✓
Lubricant	✓	✓





## CURRICULUM REFERENCES:

### UK

#### MEDICAL UNDERGRAD

*Outcomes for Graduates 2018, General Medical Council, p.16*

#### PRE-REG NURSE

*NMC Future nurse: Standards of proficiency for registered nurses, 17 May 2018, p.35*

#### POST GRAD NURSE

*RCN Competencies - Advanced Nurse Practitioners, 2012, p.4*

#### PHYSICIAN ASSOCIATE

*NHS National Practitioner Program: Matrix Specification for the Physician Associate (2016). p. 73*

#### EARLY POSTGRAD MEDICINE (FOUNDATION 1/2)

*The UK Foundation Programme Curriculum 2016. p.9*

#### CORE MEDICAL TRAINING

*Specialty Training Curriculum for Core Medical Training, 2013 p.86*

#### GP TRAINING

*Guidance for Satisfactory Progression at ARCP panels for GP Specialty Trainees, 2016 p. 6*

#### GENERAL SURGERY SPECIALIST TRAINING

*Intercollegiate Surgical Skills Curriculum, General Surgery, 2013 p.32 & p.187*

#### UROLOGY SPECIALIST TRAINING

*The Intercollegiate Surgical Curriculum, Urology Surgery, Aug 2015 p.20*

### USA

#### UNDERGRADUATE MEDICINE

*Undergraduate Medicine - AAMC (2008) - Recommendations for Preclerkship Clinical Skills Education for Undergraduate Medical Education p.26*

#### POSTGRADUATE NURSING

*National Organisation of Nurse Practitioner Faculties, Nurse Practitioner Core Competencies Content , 2017 Independent Practice Competencies p.14*

#### FAMILY MEDICINE

*AAFP Recommended Curriculum Guidelines for Family Medicine Residents, Men's Health Reprint No. 257, 2016*

#### GENERAL SURGERY SPECIALIST TRAINING

*SCORE Curriculum Outline for General Surgery (2017-28): p.19*

### AUSTRALIA & NEW ZEALAND

#### UNDERGRADUATE MEDICAL

*Medical Deans' Clinical Assessment Blueprints for the medical graduate*

#### POSTGRADUATE MEDICAL – JUNIOR DOCTOR

*CPMEC Australian Curriculum Framework for Junior Doctors v 3.1, 2012, p.7*

#### GP TRAINING

*RACGP Curriculum for Australian General Practice 2016, CS16 Core Skills unit, p.23*



## FOR DEMONSTRATION OR TO ORDER

Limbs & Things Ltd. Sussex Street, St Philips, Bristol, BS2 0RA, UK

T: 00 44 (0)117 311 0500 E: sales@limbsandthings.com

[limbsandthings.com](http://limbsandthings.com)