

Product Recycling Information: CFPT Product Range

Below is a pictorial guide to components within the Limbs & Things CFPT product range, their material/s and recyclable status.

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Part Description	Model	Material	Recyclable	Disposal Information
Small Carry Case 	All CFPT Product No. 40210	PET, Foam		Widely recycled
CFPT Modules 	All CFPT Modules	Silicone, foam, Lycra		Specialist recycling or landfill
Comfigel 42g Tube Lubricant 	All CFPT	Tube (dependant on material)		Recyclable, check supplied material for details
		Lubricant (water soluble), cardboard		Lubricant - safe to dispose Cardboard - Widely recycled
CFPT Support Base 	All CFPT	Silicone		Specialist recycling or landfill
CFPT Pelvis (with 16mm Plugs) 	All CFPT	Silicone		Specialist recycling or landfill
CFPT Abdominal Insert 	All CFPT Product No. 60910/60936	Silicone, Foam		Specialist recycling or landfill