


Product Recycling Information: Chest Drain & Needle Decompression Trainer Product Range

Below is a pictorial guide to components within the Limbs & Things Chest Drain & Needle Decompression Trainer product range, their material/s and recyclable status.



Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Part Description	Model	Material	Recyclable	Disposal Information
Chest Drain Torso 	All Chest	Mixed	 (See component breakdown below)	Specialist recycling or landfill
Torso Panels		MDPE		Recyclable, non-standard recycling
Bottle Tops		HDPE		Recyclable, non-standard recycling
Silicone tube		Silicone		Specialist recycling or landfill
Large Polyester Tube		Polyester		Recyclable, non-standard recycling
Tee Tapper Connectors		Polypropylene		Widely recycled

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Chest Drain Back Plate 5mm ABS		ABS		Recyclable, non-standard recycling
Panel Mount Coupler Nitto S		Polyacetal Resin		Specialist recycling or landfill
Hosetail		Nylon		Specialist recycling or landfill
In-line hose barb				
A2 S/Steel M6 x 12mm Button Head Screw		Stainless Steel		Widely recycled
A2 S/Steel M6 x 16mm Socket Head				
P/Plast30 Torx Pan				
Advanced Chest Drain Pads 	Adv Chest Product No. 60232/61232	Silicone		Specialist recycling or landfill
Standard Chest Drain Pads 	Std Chest Product No. 60234/61234	Silicone, Flexocel		Specialist recycling or landfill
Chest Drain skin 	All Chest	Silicone		Specialist recycling or landfill
Needle Decompression Pump 	All Chest	Silicone, Polyacetal Resin, Nylon		Specialist recycling or landfill
		Polypropylene (Check Valve Barb)		Widely recycled

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Respiratory Swing Pump 	All Chest	Silicone, Polyacetal Resin		Specialist recycling or landfill
		Polypropylene (Reduct Connector)		Widely recycled
Needle Decompression Pads 	All Chest Product No. 60231/61231	Silicone		Specialist recycling or landfill
Chest Drain Rib Frames (Left & Right) 	All Chest Product No. 60235/60237	ABS		Recyclable, non-standard recycling
Chest Drain Ribs 	All Chest Product No. 60240	ABS		Recyclable, non-standard recycling
Chest Drain Slotted Retaining Knobs 	All Chest	Stainless Steel		Widely recycled
Chest Drain Lung 	All Chest Product No. 60238	RTV 10A 0015 & 40A		Specialist recycling or landfill
Chest Drain Diaphragm 	All Chest Product No. 60239	RTV 10A 0015 & 40A		Specialist recycling or landfill

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Slotted Nut Key with Keyring 	All Chest	Aluminium		Widely recycled
		GPPS Plastic		Recyclable, non-standard recycling
Chest Drain Augmented Reality Mat 	All Chest Product No. 60132	Digiflor Tactile fire-rated slip-rated cushioned vinyl flooring 1000gsm		Recyclable, non-standard recycling
Packaging 	All Chest	Cardboard		Widely recycled
Medium Carry Case 	All Chest Product No. 40210	PET, foam		Widely recycled

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.