

Product Recycling Information: Abdominal Examination Trainer


Product Range

Below is a pictorial guide to components within the Limbs & Things Abdominal Examination Trainer product range, their material/s and recyclable status.

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Part Description	Model	Material	Recyclable	Disposal Information
AE Skin 	All Abdo Product No. 60006/60057	Silicone		Specialist recycling or landfill
AE Organ Locator 	All Abdo Product No. 60008	Silicone		Specialist recycling or landfill
AE Abdominal Wall 	All Abdo Product No. 60009	Foam		Specialist recycling or landfill
AE Livers 	All Abdo Product No. 60010	Silicone		Specialist recycling or landfill
AE Spleens 	All Abdo Product No. 60011	Silicone		Specialist recycling or landfill
AE Kidney 	All Abdo Product No. 60012	Silicone		Specialist recycling or landfill
AE Bladder 	All Abdo Product No. 60013	Silicone		Specialist recycling or landfill

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

AE Pathologies 	All Abdo Product No. 60014	Silicone		Specialist recycling or landfill
AE Aortas 	All Abdo Product No. 60015	Silicone		Specialist recycling or landfill
Carry Case (x2) 	All Abdo	PET		Widely recycled
Augmented Reality Mat 	All Abdo Product No. 60134	Digiflor Tactile fire- rated slip-rated cushioned vinyl flooring 1000gsm		PVC, recyclable, non-standard recycling
Packaging 	All Abdo	Cardboard		Widely recycled

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Part Description	Model	Material	Recyclable	Disposal Information
Abdominal Base Unit 	Abdominal Exam 60005	Mixed		(See component breakdown below)
AE Chest		Icorene 3560 Icbag 20 P1200		Widely recycled
Polyplast30 Torx Pan 4 x 14 SS				
Polyplast30 Torx Pan 4 x 25 SS				
AE Breathing Mechanism Piston Assembly				
3x10 Torx Pan SS 202101452K3010S				
AE Guide Rod				
A2 M6 X 12mm Button Head Screw (Bolt)		Stainless Steel		Widely recycled
M4 Button Head Screws				
M4 Washer				
Washer-M4 Shakeproof				
M8 x 30mm Socket Button Head Screw				
M8 x 35mm Socket Button Screw				
AE Breathing Mechanism Core Assembly		Stainless Steel/ Delrin 150		Widely recycled

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

6297 Ashby Step Two Rubber Foot	Rubber		Specialist recycling or landfill
AE Ribs	OnFlex S EF 80A-3S1977		Specialist recycling or landfill
AE Breathing Mechanism	Mixed		Specialist recycling or landfill
AE Interface Plate	CYCLOY C6200		Specialist recycling or landfill
AE Breathing Mechanism Casing			
AE Speaker Casing			
AE Mounting Plate			
AE MP3 Plyr TDM 156	Mixed/Electrical		Specialist recycling - e-waste & battery recycling facilities
32GB SanDisk ultra Fit USB Drive			
Battery Holder - 4AA			
Jack Socket 106-874 JST Cable Assembly			
RS Pro Butt Wire Splice Connector, Yellow Insulated 26 → 22 AWG			
DC Socket 2.1mm JST Cable Assembly			
Visaton BF 32- 80hm Speaker			
Assmann AK-AV100 3.5mm Cable 1.5m			
JST-PH Jumper Assembly (2 Wire)			
SD Card Slot Cover			
5V DC Power Supply (High Altitude)			

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

AE Core and Piston Assembly	Mixed		Specialist recycling or landfill
Deep Grove Ball Bearing	Steel		Widely recycled
M3 Nylock Nut BZP			
Screw-3 X 16mm-Pozi-Pan Head-Z.P.	Steel/Zinc		Widely recycled
Panel mount coupler Nitto SPC-06PHB (NBR) Male	Polyacetal resin		Specialist recycling or landfill
AE Bulkhead Plate	Aluminium		Widely recycled
AE Organ Mounting Plate Assembly	Thermoplastics		Recyclable, non-standard recycling
Plastic Wing Knob			
AE Rib Sound Damping Spacer	Neoprene		Specialist recycling or landfill
AE Velcro Patch	Velcro		Specialist recycling or landfill
Black Snap Together Rivet - 9mm	Nylon		Specialist recycling or landfill
White Snap Together Rivet - 9mm			
Nylon Shoulder Washer 4mm x 7mm			
Topran Crank Handle			
AE Guide Rod Casing			

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

AE Spine	PU GD1080A Thermoplastic		Recyclable, non- standard recycling
AE Pelvis Assembly			
Silicone Tube 4.8mmID x 8mm	Silicone		Specialist recycling or landfill
AE Battery Cover			
AE Rectus Bag  All Abdo Product No. 60007 Mixed  (See component breakdown below)			
AE Rectus Bag			
AE Rectus Muscle	Silicone		Specialist recycling or landfill
PG10 Lycra Reinforcing Ring			
AE Rectus Muscle	Foam		Recyclable, non- standard recycling
Polybag	Polyurethane		Recyclable, non- standard recycling
AE Distension Set  All Abdo Product No. 60016 Mixed  (See component breakdown below)			
AE Distension Foam	Flexocel		Specialist recycling or landfill
AE Gaseous Distension Bag & Pump	Silicone		Specialist recycling or landfill

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.