








Product Recycling Information: Catheterisation Product Range

Below is a pictorial guide to components within the Limbs & Things Catheterisation product range, their material/s and recyclable status.




Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Part Description	Model	Material	Recyclable	Disposal Information
Augmented Reality Mat for Catheterisation 	All Cath Product No. 60131	Digiflor Tactile fire-rated slip-rated cushioned vinyl flooring 1000gsm	 3 PVC	Recyclable, non-standard recycling
Barbed Bottle Top 	All Cath	Sabic C6200 PC/ABS	 7 OTHER	Recyclable, non-standard recycling
Bladder Container 	All Cath	Wonderlite PC-111	 7 OTHER	Recyclable, non-standard recycling
Bladder Container Bung 	Male & Female Models	FST-7040		Specialist recycling or landfill
Bladder Diaphragm 	All Cath Product No. 60866	Silicone		Specialist recycling or landfill
Bladder Valve 	All Cath Product No. 60882	Silicone		Specialist recycling or landfill

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Cardboard Tube 	All Cath	Cardboard		Widely recycled
Catheterisation Lubricant Bottle 	All Cath Product No. 60892	PET Bottle		Widely recycled
Fluid Jug 	All Cath	Wonderlite PC-110		Recyclable, non-standard recycling
Fluid Jug Lid 	All Cath	TPE TC9GPN 89ShA		Recyclable, non-standard recycling
Foreskins for Catheterisation 	Male Models Product No. 60859/60878	Silicone		Specialist recycling or landfill
Kidney Dish 	All Cath	Polypropylene		Widely recycled
Metal Fluid Stand 	All Cath	S235JR Mild Steel Powder Coated		Widely recycled
Module Casing 	All Cath	Wonderlite PC-110		Recyclable, non-standard recycling
Perineum (Male & Female) 	All Cath	Silicone		Specialist recycling or landfill

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Pinch Clamp 	All Cath	Nylon		Recyclable, non-standard recycling
Pubic Bone 	All Cath	Sabic C6200 PC/ABS		Recyclable, non-standard recycling
Stand 	All Cath	Silicone		Specialist recycling or landfill
Suprapubic Bladder 	Suprapubic Models	Silicone		Specialist recycling or landfill
Suprapubic Connector 	Suprapubic Models	Acetal Resin		Specialist recycling or landfill
Suprapubic Insertion Bungs 	Suprapubic Models	Mixed – Silicone & Fabric		Specialist recycling or landfill
Suprapubic Perineum 	Suprapubic Models	Silicone		Specialist recycling or landfill
Torso (Lower) 	All Cath	Silicone		Specialist recycling or landfill

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.

Urethra Support Plate 	Male Models	Polypropylene		Widely recycled
Water Extension Tubing 	All Cath Product No. 60867	Silicone		Specialist recycling or landfill

Available recycling and disposal options will depend on local facilities and infrastructure. Check with your waste management provider for more information on suitable recycling, recovery, and disposal routes.