

Injection Trainer

Part No: 00310

**Learning Has Never
Been More Lifelike**

Epidermis layer may be peeled back, to release the subcuticular injection liquid, and replaced in seconds


Multiple tissue layers: epidermis, dermis, fat & muscle

Trainer is easy to assemble and disassemble

All components are latex free

Trainer can be strapped to a trainee's arm or thigh for role playing

Liquid injected into the pad can be squeezed out



Injection Trainer

Part No: 00310

Multi-layered soft tissue pad for the practice of intradermal, subcutaneous, and intramuscular tissue injection techniques.



Skills Gained:

- Professional-to-patient communication
- Subcutaneous injection
- Intradermal injection and Intramuscular injection
- Management of tissue

Product Qualities:



- Multiple tissue layers: epidermis, dermis, fat & muscle
- Long-lasting replaceable epidermis
- Trainer can be strapped to a trainee's arm or thigh for role playing



- Epidermis layer may be peeled back, to release the subcuticular injection liquid, and replaced in seconds
- Liquid injected into the intramuscular pad can be squeezed out
- Trainer is easy to assemble and disassemble



- All components are latex free



You'll Receive:

- 00311 1 Skin Pad & Muscle Block
- 00320 1 Epidermis for Injection Trainer (Pack of 2)
- 1 Strap-on Base Plate
- 1 Soft Tissue Tray with Liner

Packed weight (kg) & dimensions (mm) 600 g 210 l 160 w 110 h



limbsandthings.com

Limbs & Things Australia Pty Ltd, 92 Wedgewood Road, Hallam, Victoria 3803, Australia
E: info@limbsandthings.com.au T: +61 3 9708 6511 F: +61 3 9708 6566